

Product Datasheet

Phospho-Gli1 (T101/ S102/S104) Rabbit Polyclonal Antibody (Cy7) (orb1083786)

Catalog Number	orb1083786
Category	Antibodies
Description	Phospho-Gli1 (T101/ S102/S104) Rabbit Polyclonal Antibody (Cy7) is a Cy7 conjugated antibody targeting GLI1. This antibody is suitable for IF. It exhibits reactivity with Human, Mouse samples.
Target	GLI1
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	Cy7
Reactivity	Human, Mouse
Predicted Reactivity	Bovine, Canine, Equine, Porcine, Rat, Sheep
Concentration	1mg/1ml
Buffer/Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Purification	Affinity purified by Protein A
Immunogen	KLH conjugated Synthesised phosphopeptide peptide derived from human Gli1 around the phosphorylation site of T101/S102/S104 IR(p-T)(p-S)P(p-S)SL
Tested applications	IF

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Dilution range	IF=1:100-500
Antibody Type	Primary Antibody
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)