

Product Datasheet

Sphingosine 1 Phosphate 4 Receptor peptide (orb108398)

Catalog Number	orb108398
Category	Proteins
Description	Sphingosine 1 Phosphate 4 Receptor peptide
Conjugation	Unconjugated
Reactivity	Human
Form/Appearance	Provided in solution at 1 mg/ml in phosphate buffered saline.
Concentration	See vial for concentration
Tested applications	FA
Application notes	For use with EDG6 (S1P4) (CT) polyclonal antibody (orb108773).
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles. See expiration date on vial
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)