

Product Datasheet

Sphingosine 1 Phosphate 4 Receptor peptide (orb108472)

Catalog Number	orb108472
Category	Proteins
Description	Sphingosine 1 Phosphate 4 Receptor peptide
Form/Appearance	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Concentration	See vial for concentration
UniProt ID	O95977
Tested applications	Blocking
Application notes	For use with S1P4 polyclonal antibody (orb108892). Optimal concentration should be evaluated by serial dilutions.
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles. See expiration date on vial
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)