

Product Datasheet

Human IL-1 β ELISA Kit (orb1085827)

Catalog Number	orb1085827
Category	Assays and Kits
Description	The Human IL-1 β ELISA Kit features a dynamic range of 7.81- 500 pg/mL and is suitable for Cell culture supernates, serum, and plasma. This product offers a sensitivity of 4 pg/mL.
Reactivity	Human
Range	7.81- 500 pg/mL
Sensitivity	4 pg/mL
Sample Types	Cell culture supernates, serum, and plasma.
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice.).
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)