


Product Datasheet

Human IFN- γ ELISA Kit (orb1085865)

Catalog Number	orb1085865
Category	Assays and Kits
Description	The Human IFN- γ ELISA Kit features a dynamic range of 15.60 - 1,000 pg/mL and is suitable for Cell culture supernates, serum, and plasma. This product offers a sensitivity of 7 pg/mL.
Reactivity	Human
Range	15.60 - 1,000 pg/mL
Sensitivity	7 pg/mL
Sample Types	Cell culture supernates, serum, and plasma.
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice).
Note	For research use only
Expiration Date	Please enquire.



Standard curve for orb1085865 ELISA kit.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558