

Product Datasheet

Swine IL-1 β ELISA Kit (orb1086029)

Catalog Number	orb1086029
Category	Assays and Kits
Description	The Swine IL-1 β ELISA Kit features a dynamic range of 62.50 - 4,000 pg/mL and is suitable for Cell culture supernates, serum, and plasma. This product offers a sensitivity of 30 pg/mL.
Reactivity	Porcine
Range	62.50 - 4,000 pg/mL
Sensitivity	30 pg/mL
Sample Types	Cell culture supernates, serum, and plasma.
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice.).
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)