

## Product Datasheet

### Human VCAM-1 ELISA Kit (orb1086081)

<b>Catalog Number</b>	orb1086081
<b>Category</b>	Assays and Kits
<b>Description</b>	The Human VCAM-1 ELISA Kit features a dynamic range of 15.60 -1,000 pg/mL and is suitable for Cell culture supernates, serum, and plasma. This product offers a sensitivity of 7 pg/mL.
<b>Reactivity</b>	Human
<b>Range</b>	15.60 -1,000 pg/mL
<b>Sensitivity</b>	7 pg/mL
<b>Sample Types</b>	Cell culture supernates, serum, and plasma.
<b>Storage</b>	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice.).
<b>Note</b>	For research use only
<b>Expiration Date</b>	Please enquire.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)