

Product Datasheet

Human Peroxiredoxin 2 ELISA Kit (orb1086185)

Catalog Number	orb1086185
Category	Assays and Kits
Description	The Human Peroxiredoxin 2 ELISA Kit features a dynamic range of 781.25 - 50000 pg/mL and is suitable for Cell culture supernates, serum, and plasma. This product offers a sensitivity of 390 pg/mL.
Reactivity	Human
Range	781.25 - 50000 pg/mL
Sensitivity	390 pg/mL
Sample Types	Cell culture supernates, serum, and plasma.
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice.).
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)