

Product Datasheet

Human IL-18 BPa ELISA Kit (orb1086193)

| | |
|------------------------|---|
| Catalog Number | orb1086193 |
| Category | Assays and Kits |
| Description | The Human IL-18 BPa ELISA Kit features a dynamic range of 93.75-6000 pg/mL and is suitable for Cell culture supernates, serum, and plasma. This product offers a sensitivity of 47 pg/mL. |
| Reactivity | Human |
| Range | 93.75-6000 pg/mL |
| Sensitivity | 47 pg/mL |
| Sample Types | Cell culture supernates, serum, and plasma. |
| Storage | Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice.). |
| Note | For research use only |
| Expiration Date | Please enquire. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)