

## Product Datasheet

### Mouse CXCL11/I-TAC ELISA Kit (orb1086279)

<b>Catalog Number</b>	orb1086279
<b>Category</b>	Assays and Kits
<b>Description</b>	The Mouse CXCL11/I-TAC ELISA Kit features a dynamic range of 62.50 - 4,000 pg/mL and is suitable for Cell culture supernates, serum, and plasma. This product offers a sensitivity of 31 pg/mL.
<b>Reactivity</b>	Mouse
<b>Range</b>	62.50 - 4,000 pg/mL
<b>Sensitivity</b>	31 pg/mL
<b>Sample Types</b>	Cell culture supernates, serum, and plasma.
<b>Storage</b>	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice.).
<b>Note</b>	For research use only
<b>Expiration Date</b>	Please enquire.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)