

Product Datasheet

Mouse Cytokeratin 7, clone OV/TL 12/30 (monoclonal) Antibody (orb1087758)

Catalog Number	orb1087758
Category	Antibodies
Description	Mouse anti-Cytokeratin 7, clone OV/TL 12/30 (monoclonal)
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG1, kappa
Conjugation	Unconjugated
Reactivity	Human
Concentration	Greater than or equal to 67 mg/L
Buffer/Preservatives	Tissue culture supernatant with Sodium azide
Immunogen	Please contact us for more information.
Tested applications	IHC-P
Dilution range	IHC-P: 1:50
Clone Number	OV-TL 12/30
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Expiration Date

12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)