

Product Datasheet

Mouse Cytokeratin 7, clone OV-TL 12/30 Antibody (orb1087974)

Catalog Number	orb1087974
Category	Antibodies
Description	Mouse anti-Cytokeratin 7, clone OV-TL 12/30
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG1-k
Conjugation	Unconjugated
Reactivity	Human
Buffer/Preservatives	Tris Buffer, pH 7.3-7.7, containing 1% rAlbumin and 0.1% Sodium Azide
Immunogen	Please contact us for more information.
Tested applications	IHC-P
Clone Number	OV-TL 12/30
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)