

Product Datasheet

Human IL-12(p70) ELISA Kit (DIY Antibody Pairs) (orb1098132)

Catalog Number	orb1098132
Category	Assays and Kits
Description	Human IL-12 (p70) ELISA Kit (DIY Antibody Pairs) (antibody pairs and standards for assay development). Quantitate Human IL-12 (p70) in cell culture supernatants and serum.
Target	Interleukin-12 subunit alpha
Reactivity	Human
Immunogen	Expression system for standard: NS0; Immunogen sequence: (I23-S328) + (R23-S219)
UniProt ID	P29459
MW	28293 Da
Tested applications	ELISA
Specificity	Natural and recombinant human IL-12(p70)
Cross Reactivity	There is no detectable cross-reactivity with other relevant proteins.
Range	7.8 pg/ml - 500 pg/ml
Sample Types	cell culture supernatants and serum.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Expiration Date 6 months from date of receipt.



Human IL-12(p70) EZ Set ELISA Kit standard curve

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)