

Product Datasheet

Human CCL2/MCP1 ELISA Kit (DIY Antibody Pairs) (orb1098147)

Catalog Number	orb1098147
Category	Assays and Kits
Description	Human CCL2/MCP1 ELISA Kit (DIY Antibody Pairs) (antibody pairs and standards for assay development). Quantitate Human CCL2 in cell culture supernatants, serum, plasma (heparin, EDTA, citrate) and urine.
Target	C-C motif chemokine 2
Reactivity	Human
Immunogen	Expression system for standard: E.coli; Immunogen sequence: Q24-T99
UniProt ID	P13500
Tested applications	ELISA
Specificity	Natural and recombinant human MCP-1
Cross Reactivity	There is cross-reactivity with human Eotaxin, MCP-31%.
Range	15.6 pg/ml - 1,000 pg/ml
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA, citrate) and urine.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.


7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558



Human CCL2/MCP1 EZ Set ELISA Kit standard curve

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)