

Product Datasheet

Human IGFBP3 ELISA Kit (DIY Antibody Pairs) (orb1098166)

Catalog Number	orb1098166
Category	Assays and Kits
Description	Human IGFBP3 ELISA Kit (DIY Antibody Pairs) (antibody pairs and standards for assay development). Quantitate Human IGFBP3 in cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
Target	Insulin-like growth factor-binding protein 3
Reactivity	Human
Immunogen	Expression system for standard: NS0; Immunogen sequence: A28Q-291
UniProt ID	P17936
Tested applications	ELISA
Specificity	Natural and recombinant human IGFBP3
Cross Reactivity	There is no detectable cross-reactivity with other relevant proteins.
Range	156 pg/ml - 10,000 pg/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com


Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



Human IGFBP3 EZ-Set ELISA Kit standard curve

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)