

Product Datasheet

Human IGFBP2 ELISA Kit (DIY Antibody Pairs) (orb1098169)

Catalog Number	orb1098169
Category	Assays and Kits
Description	Human IGFBP2 ELISA Kit (DIY Antibody Pairs) (antibody pairs and standards for assay development). Quantitate Human IGFBP2 in cell culture supernatants, serum, plasma (heparin, EDTA) and urine.
Target	Insulin-like growth factor-binding protein 2
Reactivity	Human
Immunogen	Expression system for standard: NS0; Immunogen sequence: E40-Q328
UniProt ID	P18065
Tested applications	ELISA
Specificity	Natural and recombinant human IGFBP2
Cross Reactivity	There is no detectable cross-reactivity with other relevant proteins.
Range	62.5 pg/ml - 4,000 pg/ml
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and urine.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.


7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



Human IGFBP2 EZ-Set ELISA Kit standard curve

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)