

Product Datasheet

Mouse VEGF-B ELISA Kit (DIY Antibody Pairs) (orb1098192)

Catalog Number	orb1098192
Category	Assays and Kits
Description	Mouse VEGF-B ELISA Kit (DIY Antibody Pairs) (antibody pairs and standards for assay development). Quantitate Mouse Vegfb in cell culture supernatants, serum and plasma (heparin).
Target	Vascular endothelial growth factor B
Reactivity	Mouse
Immunogen	Expression system for standard: sf21; Immunogen sequence: Q20-A207 & Q20-R148
UniProt ID	P49766
Tested applications	ELISA
Specificity	Natural and recombinant mouse VEGF-B
Cross Reactivity	There is no detectable cross-reactivity with other relevant proteins.
Range	31.2 pg/ml - 2,000 pg/ml
Sample Types	cell culture supernatants, serum and plasma (heparin).
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com


Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Mouse VEGF-B ELISA Kit

Mouse VEGF-B EZ-Set ELISA Kit standard curve

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)