

## Product Datasheet

### Human TNFRSF17/BCMA ELISA Kit (DIY Antibody Pairs) (orb1098218)

|                            |                                                                                                                                                                                                                  |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>      | orb1098218                                                                                                                                                                                                       |
| <b>Category</b>            | Assays and Kits                                                                                                                                                                                                  |
| <b>Description</b>         | Human TNFRSF17/BCMA ELISA Kit (DIY Antibody Pairs) (antibody pairs and standards for assay development). Quantitate Human TNFRSF17 in cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA). |
| <b>Target</b>              | Tumor necrosis factor receptor superfamily member 17                                                                                                                                                             |
| <b>Reactivity</b>          | Human                                                                                                                                                                                                            |
| <b>Immunogen</b>           | Expression system for standard: NS0; Immunogen sequence: M1-A54                                                                                                                                                  |
| <b>UniProt ID</b>          | <b>Q02223</b>                                                                                                                                                                                                    |
| <b>Tested applications</b> | ELISA                                                                                                                                                                                                            |
| <b>Specificity</b>         | Natural and recombinant human TNFRSF17/BCMA                                                                                                                                                                      |
| <b>Cross Reactivity</b>    | There is no detectable cross-reactivity with other relevant proteins.                                                                                                                                            |
| <b>Range</b>               | 31.2 pg/ml - 2,000 pg/ml                                                                                                                                                                                         |
| <b>Sample Types</b>        | cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).                                                                                                                                       |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.                                                                          |
| <b>Note</b>                | For research use only                                                                                                                                                                                            |
| <b>Expiration Date</b>     | 6 months from date of receipt.                                                                                                                                                                                   |

#### Biorbyt Ltd.


7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

**Human TNFRSF17 ELISA Kit**

Human TNFRSF17/BCMA EZ-Set ELISA Kit standard curve

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)