

Product Datasheet

Mouse CXCL3 ELISA Kit (DIY Antibody Pairs) (orb1098228)

Catalog Number	orb1098228
Category	Assays and Kits
Description	Mouse CXCL3 ELISA Kit (DIY Antibody Pairs) (antibody pairs and standards for assay development). Quantitate Mouse CXCL3 in cell culture supernatants, serum and plasma (heparin, EDTA).
Target	C-X-C motif chemokine 3
Reactivity	Mouse
Immunogen	Expression system for standard: E.coli; Immunogen sequence: S32-S100
UniProt ID	Q6W5C0
Tested applications	ELISA
Specificity	Natural and recombinant mouse CXCL3
Cross Reactivity	There is no detectable cross-reactivity with other relevant proteins.
Range	15.6 pg/ml - 1,000 pg/ml
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.


7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558



Mouse CXCL3 EZ Set ELISA Kit standard curve

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)