


Product Datasheet

Mouse ROCK1 (Rhoassociated, coiled-coil containing) ELISA Kit (orb1100867)

Description	Mouse ROCK1 (Rhoassociated, coiled-coil containing) ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/ml
Target	ROCK1
Storage	2 - 8 °C for 12 months
Note	For research use only
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Assay Time	4 hours
Uniprot ID	P70335
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Mouse ROCK1 (Rhoassociated, coiled-coil containing) ELISA Kit

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+441223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)