


Product Datasheet

Human PDCD6IP (Programmed cell death 6-interacting protein) ELISA Kit (orb1101055)

Catalog Number	orb1101055
Description	Human PDCD6IP (Programmed cell death 6-interacting protein) ELISA Kit
Reactivity	Human
Range	78.125-5000pg/ml
Target	PDCD6IP
Storage	2 - 8 °C for 12 months
Note	For research use only
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Assay Time	4 hours
Uniprot ID	Q8WUM4
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Human PDCD6IP (Programmed cell death 6-interacting protein) ELISA Kit

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)