

Product Datasheet

Anti-M112-113/E1 (MCMV) antibody - clone M112-113.02 (orb1101108)


Catalog Number	orb1101108
Description	Mouse monoclonal antibody to M112-113/E1 (MCMV)
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	The immunogen consisted of MCMV virions
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Long term -20 °C, short term +4 °C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Clone Number	M112-113.02
Purity	Affinity purified
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)



WB analysis of cells: a=uninfected cells b=MCMV infected cells
using Anti-M112-113/E1 (MCMV) antibody

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)