

Product Datasheet

Anti-Nestin (orb1122137)

Description	Mouse monoclonal antibody to Nestin
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Nestin
Preservatives	Contains 0.02% Sodium Azide
Form/Appearance	Liquid
Concentration	0.5 ml
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	WB
Isotype	IgM
Clonality	Monoclonal
Antibody Type	Primary Antibody
Purity	Cell Culture
Dilution Range	End user must optimize for experimental assay
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)