

Product Datasheet

Anti-USP14 (orb1124401)

Catalog Number	orb1124401
Category	Antibodies
Description	Rabbit polyclonal antibody to USP14
Clonality	Polyclonal
Species/Host	Rabbit
Conjugation	Unconjugated
Reactivity	Human
Form/Appearance	Liquid
Concentration	100ug/100ul
Buffer/Preservatives	PBS, pH 7.4 with 0.02% sodium azide.
Purity	Affinity purified
Purification	Affinity purified
Immunogen	Rabbit polyclonal USP14(1) antibody was raised against a peptide from the C terminal residues of human USP14 protein.
Tested applications	ELISA, IHC, WB
Antibody Type	Primary Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)