

Product Datasheet

Anti-BACE1 (orb1133168)

Catalog Number	orb1133168
Category	Antibodies
Description	Rabbit polyclonal antibody to BACE1
Clonality	Polyclonal
Species/Host	Rabbit
Conjugation	Unconjugated
Reactivity	Human
Form/Appearance	Liquid
Concentration	100ug/100ul
Buffer/Preservatives	PBS, pH 7.4 with 0.02% Sodium Azide
Purity	Affinity purified
Purification	Affinity purified
Immunogen	Affinity purified
Tested applications	ELISA, IHC, WB
Dilution range	ELISA titer: 1:20,000-1:80,000 Western blot analysis: 1:500-1:1,000 Immunohistochemistry: 1:100-500 Optimal dilutions/concentrations should be determined by the end use
Application notes	ELISA titer: 1:20,000-1:80,000 Western blot analysis: 1:500-1:1,000 Immunohistochemistry: 1:100-500 Optimal dilutions/concentrations should be determined by the end use

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Specificity	Human
Antibody Type	Primary Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)