

Product Datasheet

Anti-PMS2 Antibody (orb1140221)

Catalog Number	orb1140221
Category	Antibodies
Description	Rabbit monoclonal antibody to Recombinant PMS2
Clonality	Monoclonal
Species/Host	Rabbit
Isotype	IgG, kappa
Conjugation	Unconjugated
Reactivity	Human
Form/Appearance	Tris Buffer, pH 7.3 - 7.7, with 1% rAlbumin and 0.1% Sodium Azide
Buffer/Preservatives	Tris Buffer, pH 7.3 - 7.7, with 1% rAlbumin and 0.1% Sodium Azide
Purity	Purified
Purification	Purified
Immunogen	Recombinant Protein
Tested applications	ELISA, IHC
Dilution range	IHC: 1:1000-1:3000; ELISA: 1:1000-1:2000 (Optimal dilutions or concentrations should be determined by the end user)
Application notes	1:50-1:200
Specificity	PMS2

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Antibody Type	Primary Antibody
Clone Number	IHC422
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)