

Product Datasheet

ATI 2341 (orb1140485)

Catalog Number	orb1140485
Category	Proteins
Description	CXCR4 allosteric agonist; Peptides.
Form/Appearance	Freeze dried solid
Purity	> 95% by hplc
Protein Sequence	Pal-Met-Gly-Tyr-Gln-Lys-Lys-Leu-Arg-Ser-Met-Thr-Asp-Lys-Tyr-Arg-Leu-OH
MW	2256.82 Da
Structure	Pal-Met-Gly-Tyr-Gln-Lys-Lys-Leu-Arg-Ser-Met-Thr-Asp-Lys-Tyr-Arg-Leu-OH
Modifications	N-terminal palmitoyl.
Target Areas	Chemokine receptors
Solubility (25°C)	Soluble to 1 mg/ml in water
CAS Number	1337878-62-2
Formula	C ₁₀₄ H ₁₇₈ N ₂₆ O ₂₅ S ₂
Storage	Store dry, frozen and in the dark
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)