

Product Datasheet

Recombinant Human anti-COVID-19/SARS-CoV-2 Nucleocapsid Protein Monoclonal Antibody (orb1141368)

Description	Recombinant Human anti-COVID-19/SARS-CoV-2 Nucleocapsid Protein Monoclonal Antibody
Species/Host	Human
Endotoxins	< 0.1EU/ug determined by LAL method.
Conjugation	Unconjugated
Tested Applications	ELISA, GICA, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Recombinant
Clone Number	3A8
Antibody Type	Recombinant Antibody
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)