

Product Datasheet

Human Cathepsin A ELISA Kit (orb1145718)

| | |
|----------------------------|---|
| Catalog Number | orb1145718 |
| Category | Assays and Kits |
| Description | Human Cathepsin A ELISA Kit |
| Target | Lysosomal protective protein |
| Reactivity | Human |
| Immunogen | Expression system for standard: NS0; Immunogen sequence: A29-Y480 |
| UniProt ID | P10619 |
| Tested applications | ELISA |
| Specificity | Natural and recombinant Human CTSA |
| Cross Reactivity | There is no detectable cross-reactivity with other relevant proteins. |
| Range | 312 pg/ml - 20,000 pg/ml |
| Sensitivity | 10 pg/ml |
| Sample Types | cell culture supernatants, serum and plasma (heparin, EDTA) |
| Storage | Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.) |
| Note | For research use only |
| Expiration Date | 6 months from date of receipt. |

Biorbyt Ltd.


7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Human Cathepsin A ELISA Kit

Human Cathepsin A PicoKine ELISA Kit standard curve

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)