

Product Datasheet

MTAP Rabbit Monoclonal Antibody (orb1147419)

Catalog Number	orb1147419
Category	Antibodies
Description	MTAP Rabbit Monoclonal Antibody
Target	T-cell surface glycoprotein CD3 zeta chain
Clonality	Monoclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	Unconjugated
Reactivity	Human, Mouse, Rat
Form/Appearance	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Immunogen	A synthesized peptide derived from human MTAP
UniProt ID	Q13126
MW	31 kDa
Tested applications	FC, ICC, IF, IHC, IP, WB
Application notes	WB 1:500-2000IHC 1:50-200ICC/IF 1:50-200IP 1:50FC 1:50. Restore with deionized water (or equivalent) for reconstitution volume of 1.0 mL

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom


Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Antibody Type	Primary Antibody
Clone Number	28M98
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only



Western blot analysis of MTAP expression in (1) 293T cell lysate; (2) NIH/3T3 cell lysate.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)