

Product Datasheet

HRP Conjugated Anti- β -Actin Mouse Monoclonal Antibody (1C7) (orb1147865)

Description	Featured HRP Conjugated Anti- β -Actin Mouse Monoclonal Antibody (1C7), specially designed for your Western blot as internal control.
Reactivity	Canine, Gallus, Hamster, Human, Insect, Monkey, Mouse, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Preservatives	Liquid in PBS, pH 7.4 and 50% Glycerol.
Form/Appearance	Liquid solution
Storage	Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.
Note	For research use only
Application notes	Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB (1:1000-1:3000).
Isotype	Mouse IgG1
Clonality	Monoclonal
Hazard Information	The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)