

Product Datasheet

Anti-Rat IgG H&L Chain Antibody, HRP conjugate (orb1148127)

Catalog Number	orb1148127
Category	Antibodies
Description	Anti-Rat IgG H&L Chain Antibody, HRP conjugate
Target	Rat IgG H&L Chain, HRP conjugate
Clonality	Polyclonal
Species/Host	Goat
Conjugation	HRP
Reactivity	Rat
Form/Appearance	Liquid
Concentration	1.0 mg/mL
Buffer/Preservatives	0.01% Pro-Clean 400. PBS, 0.2% rAlbumin, 0.01% Pro-Clean 400.
Purity	Purified
Immunogen	Rat IgG Host Species: Goat Specificity: Antibody is affinity-purified on rat IgG and conjugated horseradish peroxidase (HRP). Molar enzyme:antibody ratio = 4:1. This antibody may cross-react with IgG of other species.
Tested applications	ELISA, ICC, IHC, WB
Storage	This antibody is stable for at least 1 year at 2-8CC.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)