

Product Datasheet

Anti-Mouse IgM Antibody, HRP conjugate (orb1148241)

Catalog Number	orb1148241
Category	Antibodies
Description	Anti-Mouse IgM Antibody, HRP conjugate
Target	Mouse IgM, HRP conjugate
Clonality	Polyclonal
Species/Host	Goat
Conjugation	HRP
Reactivity	Mouse
Form/Appearance	Liquid
Concentration	1.0 mg/mL
Buffer/Preservatives	0.1% Pro-Clean 400. PBS, 0.2% rAlbumin, 0.1% Pro-Clean 400, conjugated to horseradish peroxidase. Molar enzyme:antibody ratio = 4:1.
Purity	Purified
Immunogen	Mouse IgM Host Species: Goat Specificity: Antibody is affinity-purified on mouse IgM. By immunoblotting and ELISA this antibody reacts specifically with mouse IgM. Cross-reactivity with mouse IgG subclasses and IgA is negligible. This antibody may cross-react with IgM of other species.
Tested applications	ELISA, ICC, IHC, WB
Storage	This antibody is stable for at least 1 year at 2-8C.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Note

For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)