

Product Datasheet

Human TACE ELISA Kit (orb1149456)

| | |
|----------------------------|---|
| Catalog Number | orb1149456 |
| Category | Assays and Kits |
| Description | Human TACE/ADAM17 (TNF alpha Converting Enzyme) ELISA Kit |
| Target | TACE/ADAM17 |
| Reactivity | Human |
| Tested applications | ELISA |
| Application notes | Sandwich |
| Range | 0.25ng/ml-8ng/ml |
| Sensitivity | 0.1ng/ml |
| Sample Types | Cells, serum, tissue, saliva, lavage fluid, urine, etc. |
| Storage | 2-8°C |
| Note | For research use only |
| Expiration Date | 6 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)