

Product Datasheet

Human ANGPTL8 ELISA Kit (orb1149461)

Catalog Number	orb1149461
Category	Assays and Kits
Description	Human ANGPTL8 (Angiopoietin Like Protein 8) ELISA Kit
Target	ANGPTL8
Reactivity	Human
Tested applications	ELISA
Application notes	Sandwich
Range	1.25ng/ml-40ng/ml
Sensitivity	0.1ng/ml
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Storage	2-8°C
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)