

Product Datasheet

Human IL-12 p40 ELISA Kit (orb1149671)

Catalog Number	orb1149671
Category	Assays and Kits
Description	Human IL-12 p40 (Interleukin 12 p40) ELISA Kit
Target	IL-12 p40
Reactivity	Human
Tested applications	ELISA
Application notes	Sandwich
Range	1pg/ml-32pg/ml
Sensitivity	0.1pg/ml
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Storage	2-8°C
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)