

Product Datasheet

Mouse IL-27 p28 ELISA Kit (DIY Antibody Pairs) (orb1152409)

Catalog Number	orb1152409
Category	Assays and Kits
Description	Mouse IL-27 p28 ELISA Kit (DIY Antibody Pairs) (antibody pairs and standards for assay development). Quantitate Mouse Il27 in cell culture supernatants, serum and plasma (heparin, EDTA).
Target	Interleukin-27 subunit alpha
Reactivity	Mouse
Immunogen	Expression system for standard: NS0; Immunogen sequence: F29-S234
UniProt ID	Q8K3I6
Tested applications	ELISA
Specificity	Natural and recombinant mouse Il27
Cross Reactivity	There is no detectable cross-reactivity with other relevant proteins.
Range	15.6 pg/ml - 1,000 pg/ml
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA)
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.


7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



Mouse IL-27 p28 EZ Set ELISA Kit standard curve

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)