

Product Datasheet

Mouse IL-1 Beta/IL-1F2/IL1B ELISA Kit (DIY Antibody Pairs) (orb1152417)

Catalog Number	orb1152417
Category	Assays and Kits
Description	Mouse IL-1 Beta/IL-1F2/IL1B ELISA Kit (DIY Antibody Pairs) (antibody pairs and standards for assay development). Quantitate Mouse Il1b in cell culture supernatants, serum and plasma (heparin, EDTA).
Target	Interleukin-1 beta
Reactivity	Mouse
Immunogen	Expression system for standard: E.coli; Immunogen sequence: V118-S269
UniProt ID	P10749
Tested applications	ELISA
Specificity	Natural and recombinant mouse IL-1 beta
Cross Reactivity	There is no detectable cross-reactivity with other relevant proteins.
Range	12.5 pg/ml - 800 pg/ml
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA)
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.


7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558



Mouse IL-1 Beta EZ Set ELISA Kit standard curve

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)