

Product Datasheet

Mouse Monoclonal human EphA1 Antibody (orb1153491)

Catalog Number	orb1153491
Category	Antibodies
Description	Mouse Monoclonal anti-human EphA1
Clonality	Monoclonal
Isotype	Mouse IgG2
Conjugation	Unconjugated
Reactivity	Human
Purification	The IgG fraction of culture supernatant was purified by Protein G affinity chromatography and lyophilized from a 0.2 um filtered solution in phosphate-buffered saline (PBS). Lyophilized samples are stable for 2 years from the date of receipt when stored at -70C
Immunogen	Human recombinant EphA1 EC-domain
Tested applications	IHC, WB
Dilution range	WB (1 to 500-1000), IHC (1 to 20-100)
Specificity	This antibody was selected for its ability to detect human EphA1.
Clone Number	MM0261-8D36
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Expiration Date 12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)