

Product Datasheet

Cyclin D3 Recombinant Mouse Monoclonal Antibody (orb1172536)

Catalog Number	orb1172536
Category	Antibodies
Description	Cyclin D3 Recombinant Mouse Monoclonal Antibody
Target	CCND3
Clonality	Recombinant
Species/Host	Mouse
Isotype	IgG1, Kappa
Conjugation	Unconjugated
Reactivity	Human, Rat
Predicted Reactivity	Human, Rat
Form/Appearance	Liquid
Concentration	1mg/ml
Buffer/Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Purification	Affinity purified by Protein G
Tested applications	ICC, WB
Dilution range	WB=1:500-1000, ICC/IF=1:100-500
Antibody Type	Primary Antibody

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States


Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Clone Number	G12E5
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.



Blocking buffer: 5% NFDM/TBST, Primary ab Dilution: 1:1000, Primary ab incubation condition: room temperature 2h, Secondary ab: Goat Anti-Mouse IgG H&L (HRP), Lysate: BRL, Protein loading quantity: 20 µg, Exposure time: 10 s, Predicted MW: 33 kDa, Observed MW: 33 kDa.



Blocking buffer: 5% NFDM/TBST, Primary ab Dilution: 1:1000, Primary ab incubation condition: room temperature 2h, Secondary ab: Goat Anti-Mouse IgG H&L (HRP), Lysate: HeLa, Daudi, Protein loading quantity: 20 µg, Exposure time: 10 s, Predicted MW: 33 kDa, Observed MW: 33 kDa.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com


Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558



Cell line: HeLa, Fixative: 4% Paraformaldehyde,
Permeabilization: 0.1% TritonX-100, Primary ab Dilution: 1:50,
Primary incubation condition: 4°C overnight, Secondary ab:
Goat Anti-Mouse IgG, Nuclear counter stain: DAPI (Blue),
Comment: Color green is the positive signal for orb1172536.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)