

## Product Datasheet

### HRP, Rabbit Anti-Goat IgG (orb1173274)

<b>Description</b>	HRP conjugated Rabbit Anti-Goat IgG (H+L) secondary antibody
<b>Reactivity</b>	Goat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IHC-P, WB
<b>Immunogen</b>	Goat IgG whole molecule
<b>Preservatives</b>	Liquid in PBS, pH 7.4, containing 1% rAlbumin as stabilizer and 50% Glycerol.
<b>Form/Appearance</b>	Liquid solution
<b>Storage</b>	Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.
<b>Note</b>	For research use only
<b>Application notes</b>	Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB (1:10000-1:100000), IHC-P (1:500-1:2000), ELISA (1:5000-1:100000).
<b>Isotype</b>	Rabbit IgG
<b>Clonality</b>	Polyclonal
<b>Hazard Information</b>	The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.
<b>Expiration Date</b>	12 months from date of receipt.

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)