

## Product Datasheet

### IMPA1 Rabbit Polyclonal Antibody (orb1174425)

|                             |                                                                      |
|-----------------------------|----------------------------------------------------------------------|
| <b>Catalog Number</b>       | orb1174425                                                           |
| <b>Category</b>             | Antibodies                                                           |
| <b>Description</b>          | Rabbit polyclonal antibody to IMPA1                                  |
| <b>Target</b>               | IMPA1                                                                |
| <b>Clonality</b>            | Polyclonal                                                           |
| <b>Species/Host</b>         | Rabbit                                                               |
| <b>Isotype</b>              | IgG                                                                  |
| <b>Conjugation</b>          | Unconjugated                                                         |
| <b>Reactivity</b>           | Human                                                                |
| <b>Predicted Reactivity</b> | Mouse, Rat                                                           |
| <b>Form/Appearance</b>      | Lyophilized power with 1% rAlbumin and 0.02% NaN <sub>3</sub> .      |
| <b>Purification</b>         | Purified by antigen affinity column.                                 |
| <b>Immunogen</b>            | Recombinant fragment corresponding to 36-163aa of human IMPA1.       |
| <b>MW</b>                   | 30KD                                                                 |
| <b>Tested applications</b>  | ELISA, IHC, WB                                                       |
| <b>Dilution range</b>       | WB: 1:200-1:1,000                                                    |
| <b>Application notes</b>    | Reconstitute in 100μL of sterile H <sub>2</sub> O with 50% glycerol. |

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

**Storage**

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

**Note**

For research use only

---

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)