

Product Datasheet

Mouse Angiotensin I Converting Enzyme 2 ACE2 ELISA Kit (orb1176641)

Catalog Number	orb1176641
Category	Assays and Kits
Description	Mouse Angiotensin I Converting Enzyme 2 ACE2 ELISA Kit
Target	Mouse Angiotensin I Converting Enzyme 2 ACE2
Tested applications	ELISA
Range	62.5--4000 pg/mL
Sensitivity	10.0 pg/mL
Storage	Store at 2 - 8°C. Keep reconstituted standard and detection Ab at -20 °C
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)