

## Product Datasheet

# Human Heat Shock Protein 47, HSP47 ELISA Kit (orb1176718)

<b>Catalog Number</b>	orb1176718
<b>Category</b>	Assays and Kits
<b>Description</b>	Human Heat Shock Protein 47, HSP47 ELISA Kit
<b>Target</b>	Human Heat Shock Protein 47, HSP47
<b>Tested applications</b>	ELISA
<b>Specificity</b>	No cross-reaction with other related substances detected
<b>Sensitivity</b>	1.5 ng/ml
<b>Sample Types</b>	50 µL final volume, dilution factor varies on samples
<b>Storage</b>	Store at 2 - 8°C. Keep reconstituted standard and detection Ab at -20 °C
<b>Note</b>	For research use only
<b>Expiration Date</b>	6 months from date of receipt.

### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)