

Product Datasheet

Goat anti-Rat IgG-iFluor 488 (orb1178827)

Catalog Number	orb1178827
Description	Goat anti-Rat IgG-iFluor 488 is conjugated with high-activity fluorescence dye.
Species/Host	Goat
Reactivity	Rat
Conjugation	iFluor488
Tested Applications	FACS, IF
Immunogen	Purified rat IgG
Preservatives	Phosphate Buffered Saline containing 1% rAlbumin, 0.01% thimerosal, pH 7.4.
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store at 4°C. Protect from light.
Note	For research use only
Clonality	Polyclonal
Purity	Immunogen affinity chromatography.
Dilution Range	IF 1:200-1:1000 FACS 1:200-1:1000
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)