

Product Datasheet

Goat anti-Rat IgG(H+L) (orb1178841)

Catalog Number	orb1178841
Category	Antibodies
Description	Goat anti-Rat IgG(H+L) is purified from goat antiserum by immunoaffinity chromatography. It reacts with rat IgG heavy and light chains.
Clonality	Polyclonal
Species/Host	Goat
Conjugation	Unconjugated
Reactivity	Rat
Form/Appearance	Liquid
Concentration	5 mg/mL
Buffer/Preservatives	Lot-dependent, please refer to product label.
Purity	Immunogen affinity chromatography.
Immunogen	Purified rat IgG
Tested applications	WB
Dilution range	WB 1:5000-1:10000
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)