

Product Datasheet

Mouse anti-HA tag mAb (orb1178854)

Description	Mouse anti-HA tag mAb is specific to YPYDVPDYA (HA) tag, which can be applied to several immunoassays, including ELISA, WB, IP, IFA, and FACS. This antibody has been identified in several sources (bacteria, yeast, and mammalian cells) of HA tag fusion protein.
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FACS, IF, IP, WB
Immunogen	Peptide: YPYDVPDYA (HA)
Preservatives	Phosphate Buffered Saline containing 50% glycerol and 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1 kappa
Clonality	Monoclonal
Clone Number	98-106
Purity	> 95% (SDS-PAGE)
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)