

Product Datasheet

Human CD25 Mouse Monoclonal Antibody (FITC) (orb1184002)

Catalog Number	orb1184002
Category	Antibodies
Description	Human CD25 Mouse Monoclonal Antibody (FITC) is an FITC conjugated antibody targeting IL2RA. This antibody is suitable for FC, ICC, IF. It exhibits reactivity with Human samples.
Target	IL2RA
Clonality	Monoclonal
Species/Host	Mouse
Isotype	Mouse IgG2a, k
Conjugation	FITC
Reactivity	Human
Form/Appearance	Liquid
Buffer/Preservatives	0.01M PBS (pH7.4), 0.5% rAlbumin, 0.02% Proclin300.
Purification	Affinity purified by Protein G
Immunogen	KLH conjugated synthetic peptide derived from human CD25
MW	55 kDa
Tested applications	FC, ICC, IF
Dilution range	Flow-Cyt=10ul/Test, IF=1:100-500, ICC/IF=1:100-500

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Antibody Type	Primary Antibody
Clone Number	B7G1
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)