

Product Datasheet

Mouse CD25 (IL2RA) Rat Monoclonal Antibody (FITC) (orb1184114)

Catalog Number	orb1184114
Category	Antibodies
Description	Mouse CD25 (IL2RA) Rat Monoclonal Antibody (FITC) is an FITC conjugated antibody targeting IL2RA. This antibody is suitable for FC. It exhibits reactivity with Mouse samples.
Target	IL2RA
Clonality	Monoclonal
Species/Host	Rat
Isotype	Rat IgG1, λ
Conjugation	FITC
Reactivity	Mouse
Predicted Reactivity	Mouse
Form/Appearance	Liquid
Buffer/Preservatives	0.01M PBS (pH7.4), 0.5% rAlbumin, 0.02% Proclin300.
Purification	Affinity purified by Protein G
MW	55 kDa
Tested applications	FC
Dilution range	Flow-Cyt=10ul/Test

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Antibody Type	Primary Antibody
Clone Number	B6H0
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)