

Product Datasheet

Human CD107a Mouse Monoclonal Antibody (FITC) (orb1184443)

Catalog Number	orb1184443
Category	Antibodies
Description	Mouse monoclonal antibody CD107a conjugated to FITC
Target	LAMP1
Clonality	Monoclonal
Species/Host	Mouse
Isotype	Mouse IgG1, k
Conjugation	FITC
Predicted Reactivity	Human
Form/Appearance	Liquid
Concentration	Liquid
Buffer/Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Purification	affinity purified by Protein G
MW	42 kDa
Tested applications	FC
Dilution range	Flow-Cyt=10ul/Test
Antibody Type	Primary Antibody

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Clone Number	3C10
Storage	Store at 2-8°C. Protect from light. Avoid freezing.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)